




erytra
eflexis[®]



Execute excellence by
expanding the power of
your lab automation



TYPING



Technical specifications

LIST OF TECHNIQUES

Full ABO/Rh blood group typing (forward and reverse)

Rh + Kell phenotype

Antibody screen

Antibody identification

Enzyme assays

Antibody titration

Direct Coombs test

Rare antigens typing

Weak D test

ABO/Rh (D) confirmation

Crossmatch test

THROUGHPUT

Maximum throughput: 50 cards/hr (400 results/hr)¹

Full ABO/Rh blood group typing: 35 samples/hr²

Antibody screening (I, II, III): 76 samples/hr

Type & Screen: 25 samples/hr³

Throughput depends on the profiles being processed:

1 Results obtained with DG Gel® Rh Pheno + Kell

2 Results obtained with DG Gel® ABO/Rh (2D) (profile: A, B, AB, D1, D2, Ctl., N/A1, N/B)

3 Results obtained with DG Gel® ABO/Rh (2D) and antibody screening (I, II, III)

CAPACITY

2 lab configuration options:

- For peak hours:
 - 23 reagents (1 rack): 17 agitated, 5 non-agitated and 1 diluent (non-agitated reagents can be extended to 22)
 - 72 samples: 3 racks x 2 holders x 12 positions
- For more comprehensive menu:
 - 46 reagents (2 racks): 34 agitated, 10 non-agitated and 2 diluents (non-agitated reagents can be extended to 44)
 - 48 samples: 2 racks x 2 holders x 12 positions

Up to 200 DG Gel cards

Washing solutions:

- 6 L DG Fluid A solution (extendable to 12 L)*
- 6 L DG Fluid B solution (extendable to 12 L)*
- 12 L waste solution*

**With optional direct drainage to laboratory*

SIZE AND WEIGHT

Width: 110 cm (without monitor) (43.3 in)

Depth: 71 cm (27.9 in)

Height: 91 cm (35.8 in)

Weight: Approximately 173 kg (381.4 lb)

POWER SUPPLY

Voltage: 100-240 V

Frequency: 50-60 Hz

Input power: 600 W

REFERENCE	PRODUCT LIST	PRODUCT DESCRIPTION
210600	Erytra Eflexis	A fully automated, mid-sized analyzer to perform the pretransfusion compatibility tests using DG Gel technology with a smart and compact design for intuitive operations
210601	Erytra Eflexis rack for samples	Rack to load samples into the Erytra Eflexis
234179	System solution container	Container designated to hold diluted wash solution (A or B) or decontamination solution
213713	Samples holder	Holder to load samples into the analyzer
210603	Erytra Eflexis card waste container	Container for card waste

The Erytra Eflexis® analyzer is CE marked in accordance with In Vitro Diagnostic Regulation (IVDR). Product registration and availability vary by country. Ask your local Grifols representative for more information.

DIAGNOSTIC
**Support
& Services**

We deliver solutions to facilitate high-quality patient care with a dedicated and responsive team of professionals.

GRIFOLS

Distributed by:



50 Valleywood Drive
Markham, Ontario L3R 6E9
800.387.9643
800.268.1150 Français
905.470.2381 Fax
inter-medico.com

Grifols International, S.A.
Parc Empresarial Can Sant Joan
Av. de la Generalitat, 152-158
08174 Sant Cugat del Vallès
Barcelona, Spain
Tel: +34 935 710 500