

# MedicoDx™ Storage

*A Solution for Any Budget*

## SLIDE/CASSETTE STORAGE SYSTEM

- Support insert forms double layer on bottom, middle and sides
- Easy-to-assemble format – saves on storage space
- Economical
- Stackable
- Dimensions (WxHxD); 12" x 4" x 11½"



	<b>IMSCSS-10</b>	<b>Cassette Storage System (10 boxes/pkg)</b>
<i>Complete pkg:</i>	10 Boxes, 10 Dividers, 80 Trays	
	<b>IMSCSS-10-NT</b>	<b>Slide-Cassette Storage System (no Trays/10 boxes/pkg)</b>
<i>Complete pkg:</i>	10 Boxes, 10 Dividers	
	<b>IMSCS-10-S</b>	<b>Slide Storage System (10 boxes/pkg)</b>
<i>Complete pkg:</i>	10 Boxes, 10 Dividers, 40 Trays	



## CASSETTE STORAGE SYSTEM

- Easy-to-assemble format – saves on storage space
- Economical
- Available with or without divider strips
- Stores over 400 cassettes/box
- Dimensions (WxHxD); 9½" x 2" x 18"

	<b>IMCSB-50-D</b>	<b>Cassette Storage System (50 boxes/pkg)</b>
<i>Complete pkg:</i>	50 Boxes, 350 Dividers	
	<b>IMCSB-50-ND</b>	<b>Cassette Storage System (no Dividers/50 boxes/pkg)</b>
<i>Complete pkg:</i>	50 Boxes	