



MaxPlus  
EMT Cooler  
- with -  
Maxify™  
Technology



ECONOMICAL



EFFICIENT



SUSTAINABLE



USA-MADE

Product Page:



## Qualified cooler for emergency medical transport & air evac operations

Holds 2 (500 mL) whole blood units

From pre-hospital to OR and ER transfusions, trauma packs and MTP, the MaxPlus EMT cooler transports blood products safely, stores efficiently and stands up to repeated in-field and in-hospital use. All with no special handling or bench-time conditioning.

### What is Maxify™ Technology?

**MAXIFY™ proprietary materials technologies** are revolutionizing the shipping of temperature-sensitive biologics, enabling our leading-edge MAXQ no-bench-time packaging solutions to provide unmatched shipping and cost efficiencies, thermal performance, and TOTAL VALUE.

- No-Bench-Time preconditioning required
- Universal (all-season) protection
- Optimized for quick deployment
- Qualified in operational ambient environments
- Proven MaxPlus Protection - eliminates loss



Product Page:



## Maximize Savings!

MaxPlus coolers combine unprecedented thermal performance and convenience to decrease lost blood product and increase savings.


SKU #	Capacity	Dimension	Temperature	Validation
EMT10V24	(x2) 500mL units	10.875"x10.5"x9.5"	1 - 6°C	21 hrs standby, 3 in extreme transport



### System Components:

- Two PCM5 Cassettes
- One PCM0 Brick

Thermal Performance Plot: Extreme Summer Example



Distributed by:



50 Valleywood Drive  
 Markham, Ontario L3R 6E9  
 800.387.9643  
 800.268.1150 Français  
 905.470.2381 Fax  
[inter-medico.com](http://inter-medico.com)



[www.PackMaxQ.com](http://www.PackMaxQ.com)