



FDA

**Cleared for
Point of Care**



HemoScreen

**Accessible & Simplified Lab-Quality CBC
for Community Health**

Instant Diagnostics for Better Patient Care


Distributed by:



50 Valleywood Drive
Markham, Ontario L3R 6E9
800.387.9643
800.268.1150 Français
905.470.2381 Fax
inter-medico.com

HemoScreen empowers community health clinics, urgent care centers, and primary care providers with immediate access to lab-accurate CBC testing.

Its intuitive, maintenance-free design eliminates the need for external reagents or specialized staff, streamlining workflows and improving patient management without reliance on off-site laboratories.



- Portable, plug-and-play, and start up within seconds
- Lab-accurate 5-part diff and 20 CBC parameters results in 5 minutes
- NO maintenance, NO calibration, and NO external reagents required
- 3 Simple & Easy steps, using a drop of venous or capillary blood

Seamless Integration - Maximum Efficiency

Designed for point-of-care use, HemoScreen features a compact footprint and an effortless 3-step operation. With no daily maintenance or calibration required, it frees up valuable time for healthcare teams, allowing them to focus on delivering fast, high-quality care.



Reliable and easy to use, HemoScreen brings advanced diagnostics directly to the provider's office, ensuring timely results and better health outcomes.



See instructions
for use



PixCell Medical Ltd. | 6 Hayezira St., POB 113, Yokneam Illit, Israel
Office: +972-4-959-3516 | Fax: +972-4-959-3518 | Mail: info@pixcell-medical.com

Rx only